



# TALC CANTAL 45-85

SiO <sub>2</sub>	Al <sub>2</sub> O <sub>3</sub>	MgO	CaO	Fe <sub>2</sub> O <sub>3</sub>	KNa <sub>2</sub> O	MnO	L.O.I.
43.04	1.23	25.65	13.93	0.51	0.92	0.03	15.05

COMPARE WITH NYTAL 100HR (BEFORE 2009):

55.20	0.31	30.00	8.42	0.16	0.34	0.17	5.41
-------	------	-------	------	------	------	------	------

Not an exact substitute. For every 100g Nyal 100HR replaced by 100g Cantal, add 12g silica (for a total of 112 g) to maintain CaMgO: SiO<sub>2</sub> balance.